



Recombinant Schizosaccharomyces pombe Probable uroporphyrinogen-III C-methyltransferase (SPCC1739.06c)

Product Code	CSB-MP530396SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74468
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MIVSWSLKDA VVVVVGSGKN AYRRVNLCLS ENSSKIFWFC RSSADGGFDV GHLQLDAKKS DAIEILPISE FDPVRSLTTL GKEETDFLVD AVFVSESNNH EKEILHRTCK RYRIPLNIID NPSLCSFTLP ATWSEPPLQI SLSTSSNGCR LAQRLLRHVV SSLPSGMPEA IERFGRVRAI TKTPEKKQWI NHVSDFWPLE KLVRMTEDDL LAIVSDNFSL PLSSSQSSSL ANSYESLSTT LDKPSLTLDP EAFPTHKRGS IALIGSGPGS PDLLTVAARK AIMKADYVLA DKLVP EAVLQ LIPRHTPLFI ARKFPGNADK AQDELHQVAF DALSRGDYVV RLKQGDPYIY GRGGEEYLFF TQHGYVPTVI PGISSALMAP ISAGIPVTHR GVADQFLVCT GTGQKGSMPK IPSFVPTQTT VFLMALHRLE ILVQALIESG WPRVLPVCIA ERVSCPDQRF IFSTLEDVVE EYNKYESLPP GLLITGYSCN TLRNTA
Source	Mammalian cell
Target Names	SPCC1739.06c
Protein Names	Recommended name: Probable uroporphyrinogen-III C-methyltransferase EC=2.1.1.107
Expression Region	1-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.