



Recombinant *Schizosaccharomyces pombe* Putative mannan endo-1,6-alpha-mannosidase C970.02 (SPCC970.02)

Product Code	CSB-EP530406SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74556
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	I SVDLNDTSSV DLATSLVADG LLNYYAGQHK GGTIGMFLPP AYWWEAGAAW NGLLNRYIAT GNSTYNELVK TSMLYQSGED SDYMPSNYTT SEGNDQAFW GLTVISAAEA NFSNPAADep QWLELAQAVF NQQVTRWDTD HCNGGLRWQI TEFNSGYNYK NTVSNGAFFQ LAARLARFTD NDTYAEWANV AYDWSQRIGF IQEDYTVFDG SSIKDNCSI EITQWTYNIG LYMAGAAYMY NYTNSTVWKT RVEGFANKTA KTFKKDIMF EPVCEIALSC NYDQTSFKGF LTRFMVYTAQ MAPFTAPLLE PLLISTAKAA AGACCGGYDG VTTCGVQWWWN NDTWDGLYGL GEQMSALEAI QAPLLLKSLQ VFKASNGGSS TGDPNAGLYT APVSFANKNF ENLRKHWMLL GFFLLVPTLV LY
Source	E.coli
Target Names	SPCC970.02
Protein Names	Recommended name: Putative mannan endo-1,6-alpha-mannosidase C970.02 EC= 3.2.1.101 Alternative name(s): Endo-alpha-1->6-D-mannanase C970.02
Expression Region	20-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.