



Recombinant Schizosaccharomyces pombe Ubiquinone biosynthesis protein coq7 (coq7)

Product Code	CSB-YP530423SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74826
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MLSRRLIPI ATKSVSPFLR SSITCRRISM GRKELSEKDS NILDSVIRVD QAGELGANQI YKGQHFILQF TDPKVAPTIQ HMWDQEKYHL ATFDNYVLKN RVRPTFLRPF WDIAGFALGA GTALLGTKAA MACTEAVETV IGGHYNDQLR ETAHLENKAP EFKEIRSHLA EFRDDELEHL NTAVEGWNAK EAPAHALLTN AIQMGCKAAI WMCKRF
Source	Yeast
Target Names	coq7
Protein Names	Recommended name: Ubiquinone biosynthesis protein coq7 Alternative name(s): Coenzyme Q biosynthesis protein 7
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.