



Recombinant *Drosophila bakoue* Alpha-amylase-related protein (Amyrel)

Product Code	CSB-EP530462DLB-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O77019
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila bakoue</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	QHNPRWWG NR NTIVHLFEWK WSDIAEECET FLAPRGFAGV QVSPVNENII SAGRPWWERY QPISYKLTR SGNEEEFADM VRRCNVDGIR IYVDVLLNHM SGDFDGVAVG TAGTEAEPSK KSFPQVPYTA QDFHPSCEIT DWNNRFQVQE CELVGLKDLN QHSDYVRSKL IEFDLHLIEL GVAGFRVDA KHMAAVDLEY IYGSLNLNI EHGFPNARP FIFQEVIDHG HETVSREEYN ELGAVTEFRF SEEIGKAFRG NNALKWLQSW GTDWGFLNSD QALTFVDNHD NQRDQGSVLN YKSPRQYKMA TAFHLAYPYG ISRMSSFAF DDHDTPPPQD AQENIISPEF DEDGACVNGW ICEHRWRQIY AMVGFKNAVR DTELSGWWDN GDNQISFCRG NKGFLAVNNN LYDLSQELNT CLPAGEYCDV ISGSLVDGAC TGKSVTVNEH GYGYIHIGSD DFDGVLALHV NAKV
Source	E.coli
Target Names	Amyrel
Protein Names	Recommended name: Alpha-amylase-related protein EC= 3.2.1.1
Expression Region	21-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.