



Recombinant Arabidopsis thaliana Ras-related protein RABH1b (RABH1B)

Product Code	CSB-EP530500DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O80501
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAPVSALAKY KLVFLGDQSV GKTSIITRFM YDKFDNTYQA TIGIDFLSKT MYLEDRTVRL QLWDTAGQER FRSLIPSYIR DSSVAVIVYD VASRQSFLNT TKWIDEVRTE RGSDEVIVLV GNKTDLVDKR QVSIEEAEAK ARELNVMFIE TSAKAGFNIK ALFRKIAAAL PGMETLSSTK QEDMVDVNLK SSNANASLAQ QQSGGCSC
Source	E.coli
Target Names	RABH1B
Protein Names	Recommended name: Ras-related protein RABH1b Short name= AtRABH1b Alternative name(s): Ras-related protein Rab6A Short name= AtRab6A
Expression Region	1-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.