



Recombinant Arabidopsis thaliana Dehydration-responsive element-binding protein 2B (DREB2B)

Product Code	CSB-MP530546DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82133
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAVYEQTGTE QPKKRKSRAR AGGLTVADRL KKWKEYNEIV EASAVKEGEK PKRKVPAKGS KKGCMKGKGG PDNSHCSFRG VRQRIWGWV AEIREPKIGT RLWLGTFPTA EKAASAYDEA ATAMYGLAR LNFPQSVGSE FTSTSSQSEV CTVENKAVVC GDVCVKHEDT DCESNPFSQI LDVREESCGT RPDSCVGHQ DMNSSLNLYDL LLEFEQQYWG QVLQEKEKPK QEEEEIQQQQ QEQQQQQLQP DLLTVADYGW PWSNDIVNDQ TSWDPNECFD INELLDLNE PGPHQSQDQN HVNSGSYDLH PLHLEPHDGH EFNGLSSLDI
Source	Mammalian cell
Target Names	DREB2B
Protein Names	Recommended name: Dehydration-responsive element-binding protein 2B Short name= Protein DREB2B
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.