



Recombinant Arabidopsis thaliana Phospholipase A1-IIbeta (At2g31100)

Product Code	CSB-BP530554DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82274
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVGDIATRWK ELSGSSKWKD LLDPLDLDR RYILHYGDMA EVGYLAFNSD RRSKYVGDSC YTKEELFART GYLKANPFYR EVTKYIYGTS SIRLPECFII KSLSREAWNK ESNWLGYIAV ATDEGKKLLG RRGIVVAWRG TIQLYEWAND FDFPLESAVM VFPGANPNDE PRVANGWLSL YTSTDPRSFR DKTSAQEQVQ EELKRLELY KNEDVTITLT GHSLGAVMSI LSAADFLHNE WPKITPSLQH SLCVTVFAFG SPQIGDRSFK RLVESLEHLH ILRVTNVPDL IPRYPVFRFT DIGEELQINT LKSEYLKRSL NLGHFHNLEA YLHGVAGTQH NQGEFKLEIN RDIALVNKGL DALEDKYLVP GHWWVLENKG MVQSDDGTWK LNGDRSKKKQ EEEDEKEENN CKFP
Source	Baculovirus
Target Names	At2g31100
Protein Names	Recommended name: Phospholipase A1-IIbeta EC= 3.1.1.-
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.