



Recombinant Arabidopsis thaliana F-box/LRR-repeat/kelch-repeat protein At2g29770 (At2g29770)

Product Code	CSB-MP530558DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82379
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVFISETSDD GSNGGDPTKN PQEEEEENLP PIPQGIPDEL IESTVLLIRR CHYPTLSLLS KTFRRVISSS ELYKSRFILN LTDSVLYALI GFSPYNTTNL YILNCNIPRN ISLHLREIKS LPPLNHGSAV VTIGYHMYVI GGHNRLHQPT SNVSIIDLRF HTSCSLPRMQ RTRVYAAAGV IDGRIYVIGG CVKRNDHWIE VFDIENRIWS SVPHHRYCNG SSLRGEFVTV SVVMQNKIYI LDSLFGFAYE PRHGTLQSLG FETQFMFLWR DPCCVIEGLL YCIDPMCVLG HAIVVYDPNE LIWRRVKGAY ILPKFCYYQC KMANFGGKLA ILGCSNSSQR GLKDVFVFEI ELENRQSGEI WGKVDLSLAIV LRSVKSPSID LFRSVTF
Source	Mammalian cell
Target Names	At2g29770
Protein Names	Recommended name: F-box/LRR-repeat/kelch-repeat protein At2g29770
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.