



# Recombinant Arabidopsis thaliana Succinyl-CoA ligase [ADP-forming] subunit beta, mitochondrial (At2g20420)

<b>Product Code</b>	CSB-BP530566DOA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O82662
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>                     LNIH EYQGAELMGK YGVNVPKGVA ASSLEEVKKA IQDVFPNESE                      LVVKSQILAG GRGLGTFKSG LKGGVHIVKR DEAEIAGKM LGQVLVTKQT                      GPQGKVVSKV YLCEKLSLVN EMYFSIILDR KSAGPLIAC KKGTSIEDL                      AEKFPDIIK VPIDVFAGIT DEDAAKVVDG LAPKAADRKD SIEQVKKLYE                      LFRKTDCTML EINPLAETST NQLVAADAKL NFDDNAAFRQ KEVFAMRDPT                      QEDPREVAAA KVDLNYIGLD GEIGCMVNGA GLAMATMDII KLHGGTPANF                      LDVGGNASEH QVVEAFKILT SDDKVKAILV NIFGGIMKCD VIASGIVNAA                      KEVALKVPVV VRLEGTNVEQ GKRILKESGM KLITADDLDD AAEKAVKALA H                 </p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	At2g20420
<b>Protein Names</b>	Recommended name: Succinyl-CoA ligase [ADP-forming] subunit beta, mitochondrial EC= 6.2.1.5 Alternative name(s): Succinyl-CoA synthetase beta chain Short name= SCS-beta
<b>Expression Region</b>	27-421
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.