



Recombinant Arabidopsis thaliana Succinyl-CoA ligase [ADP-forming] subunit beta, mitochondrial (At2g20420)

Product Code	CSB-EP530566DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82662
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	<p> LNIH EYQGAELMGK YGVNVPKGVA ASSLEEVKKA IQDVFPNESE LVVKSQILAG GRGLGTFKSG LKGGVHIVKR DEAEIAGKM LGQVLVTKQT GPQGKVVSKV YLCEKLSLVN EMYFSIILDR KSAGPLIAC KKGTSIEDL AEKFPDIIK VPIDVFAGIT DEDAAKVVDG LAPKAADRKD SIEQVKKLYE LFRKTDCTML EINPLAETST NQLVAADAKL NFDDNAAFRQ KEVFAMRDPT QEDPREVAAA KVDLNYIGLD GEIGCMVNGA GLAMATMDII KLHGGTPANF LDVGGNASEH QVVEAFKILT SDDKVKAILV NIFGGIMKCD VIASGIVNAA KEVALKVPVV VRLEGTNVEQ GKRILKESGM KLITADDLDD AAEKAVKALA H </p>
Source	E.coli
Target Names	At2g20420
Protein Names	Recommended name: Succinyl-CoA ligase [ADP-forming] subunit beta, mitochondrial EC= 6.2.1.5 Alternative name(s): Succinyl-CoA synthetase beta chain Short name= SCS-beta
Expression Region	27-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.