



Recombinant Arabidopsis thaliana Serine/threonine-protein phosphatase PP1 isozyme 7 (TOPP7)

Product Code	CSB-MP530567DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82733
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDPGTLNSVI NRLLEAREKP GKIVQLSETE IKQLCFVSRD IFLRQPNLLE LEAPVKICGD IHGQYPDLLR LFEHGGYPPN SNYLFLGDYV DRGKQSLETI CLLLAYKIKF PENFFLLRGN HESASINRIY GFYDECKRRF SVKIWRIFTD CFNCLPVAAL IDERIFCMHG GLSPELLSLR QIRDIRRPTD IPDRGLLCDL LWSDPKDVR GWGPNDRGVS YTFGSDIVSG FLKRLDLDLI CRAHQVVEDG FEFFANKQLV TIFSAPNYCG EFDNAGAMMS VSEDLTCSFQ ILKSNDKSK FSFGSRGGAK TSFPYPKVK ILSSQNSKEY N
Source	Mammalian cell
Target Names	TOPP7
Protein Names	Recommended name: Serine/threonine-protein phosphatase PP1 isozyme 7 EC= 3.1.3.16
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.