



# Recombinant Faba bean necrotic yellows virus Para-Rep C7 (C7)

<b>Product Code</b>	CSB-BP530828FEA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O91250
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Faba bean necrotic yellows virus (isolate SV292-88) (FBNYV)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPSIRATHWC FTLNFSGSIP EINWTADVQY SIWQHERVGH DHLQGYIQMK KHSVSLKMKKE LLPGAHLEMA KAPKKAIEYC QKKESAIAGP WEYGTWISSG SHKRKLMEERF EDDPEEMKLE DPGLYRRCLS RVQMKKIRES CTWNFDLRPW QDELLKTIEQ EPDDRTIIWV YGPHGGEGKS AFAKYLTLE GWWYTAGGKA TDMLYSYS LD PTCHVCIDIP RCTREEYINY SVIEQIKNRV IINTKYEPCT IRDDGHNHVH IVFCNFLPDV TRISEDRIKI INC
<b>Source</b>	Baculovirus
<b>Target Names</b>	C7
<b>Protein Names</b>	Recommended name: Para-Rep C7 Short name= Rep7 EC= 2.7.7.- EC= 3.1.21.- EC= 3.6.1.3 Alternative name(s): ATP-dependent helicase C7 Replication-associated protein of non-essential DNA C7
<b>Expression Region</b>	1-283
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.