



Recombinant *Schizosaccharomyces pombe* Meiotic fizzy-related protein 1 (mfr1)

Product Code	CSB-MP530900SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94423
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGDRFIPIRN VSNEFNFSFQ SFKECVLSHG SNLRRKTSQT IQRQFMELLS MELFGSQASR SRAFYYGEDK RKIEKMLDT PDRKSYSLSLSP ISPQSQDMLR QPQKPKRAFP KTPYKILDAP YLKNDFYLNLDWGQSNVLA VGLASSIYLW SAASGKVVQL HDFGATNHVT SVLWTGKGTQ LAVGTDSGVI YIWDIESTKS VRSLKGHSER VAALAWNDNT LTSGGKDEVI LHHDLRAPGC CAEMMKVHEQ EICGLQWDRS LGQLASGGND NNLVWWDYRS SRPLHKFEEH TAAVKAIGWS PHQRGILASG GGTIDRCLTI HNTLTGRLQN KLDGTGSQVCN MAWSKTSNEI VTTHGFAKNQ VSLWKYPSLK NIANLTAHTN RVLVLSMSPD GQSIVTGAGD ETLRFWKLFN KKPKEESTLI R
Source	Mammalian cell
Target Names	mfr1
Protein Names	Recommended name: Meiotic fizzy-related protein 1
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.