



Recombinant *Schizosaccharomyces pombe* Uncharacterized aminotransferase C417.11c (SPCC417.11c)

Product Code	CSB-EP530915SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94492
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MRNMVSEISD LSLAKEKSEA QFIARNPKSN DFHCKACDSL PGGNTRTVLH GAPFPIFIEA GYGSKLRDVD GHEYTDLFNE LTAGIYGHNSN PVIKKALMQG FDEIGISLGG TTTCELNYAE ALKSRFLSIE KIRFCNSGTE ANITAIIAAR KFTGKRAVIA MHGGYHGGPL SFAHGISPYN MDSQDFILCE YNNSTQFKEL VNSSQDIAAV IVEAMQGAGG AIPADKEFMQ TIQLECEKND IVFILDEVMT SRLSPGGLQQ IYCLKPDLTT LGKYLGGGLP FGAFGGRADI MSCFDPRLPG SLSHSGTFNN DTLTLTAGYV GLTELYTPEA VKRLNALGDG LRKDIESYCH NTKMSITGLG SIMNIHFTES GRVNSYNDA GEVIELKDLL WMDLLKEGFW IARRGMITLS LVLTESELEA FKLTIKAWIH RRMSLIRI
Source	E.coli
Target Names	SPCC417.11c
Protein Names	Recommended name: Uncharacterized aminotransferase C417.11c EC= 2.6.-.-
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.