



Recombinant *Schizosaccharomyces pombe* Probable anthranilate synthase component 1 (trp3)

Product Code	CSB-BP530932SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94582
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKIYPLDKQV QELAEKHKAN KIPYGVIPA DMLTPSVAYL KLNQGKKYSF ILESVTQGES VSRYSFIGSP YRILMANGKT DPLARLEREL KEVKTAPVEG LPSFSGGAVG YVSYDCIKYF EPTTEMPLED TLGLPEAMFF MTDDLVAFDH AYQTVKIISH VCIQQGRPIE EAYEAAVFKI NMLKKKLESP EIPLPEQKKV HLGYEAKSNV GEDGYKAFVS NLKEHIFNGD IFQAVPSQRI ARRTDLHPFN LYRHLRTVNP SPYMFYIHCD DFDIIGASPE LLVKSEHGRI INHPIAGTVP RGKTKEEDEA YAKDLLASVK DRAEHVMLVD LARNDVSRVC DLDTTSVDKL MTIEKFSHVQ HLVSQVSGVL RPDKTRFADF RSIFPAGTVS GSPKVRAIQL VYGLEKEKRG IYAGAVGRWG YEDDNMDTCI AIRTMVYKDG TVYLLQAGGGI VFDSDEQDEY VETLNKLRNSN VTAIEETEKL YAEENSSA
Source	Baculovirus
Target Names	trp3
Protein Names	Recommended name: Probable anthranilate synthase component 1 EC=4.1.3.27 Alternative name(s): Anthranilate synthase component I
Expression Region	1-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.