



Recombinant Schizosaccharomyces pombe Uncharacterized protein C1289.06c (SPBC1289.06c)

Product Code	CSB-MP530938SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94615
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	NSQNL LKQLNDEMKA RKYTETVATF SSLKPFGLD SQTINRYILF LVDRIKMLNG RGNVDEATLD KLDDILRYAI EHQEIASARF WRIMLQSYID LNLFDKASLI ADMSLSHMEF LPKDERVLGN LYISALQAKI LGGASFEQCG KIGSAIHEQL EGKEVVNELV AVYLIYTVFK DQGISSKAHK ALVQFNGLTS FHSDVIISVF VNKGLVDQAA AYLKNSNLNE RNLPTIYTTV WLLQRLFEAH HSLDPLL TIF DYYLSVSPKD ITRLTNAI LS MKQFERDR DIKKATDFIT TFINKMSDVK EFKPSISTAN TLFSIASRLK DVKWLSAGFD MIDKYGLKPT HVTYRSLLKA YCLLPSTCEQ ISQAWVNLEY RLEAISISVA DKEINLLKDC ILSQPDRDDQ QSCLQFLNMV LSKYGKHIRS P
Source	Mammalian cell
Target Names	ppr8
Protein Names	Recommended name: Uncharacterized protein C1289.06c
Expression Region	56-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.