



Recombinant Schizosaccharomyces pombe LIN1-like protein (SPBC83.09c)

Product Code	CSB-EP530955SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94693
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MKRTLRLNPGN SSTVGDVHDV FFEYDVSNVG RKRPRTKKEG YYESESEEDE DQILNKEKKE GQSEDMFSDT SEDEKRTLPN DEAQKRRDFI ENGDAERLAH KGLRNKEVLN DDSDDDDNG KYSKLRYEDI EGQEDTNQAN DLDADEEGSE ISVPSSPKRM SFNLKEDMEE GDFDENGFI RKNYDPESQY DAWLNGSVSN KKSIAAAREA EQKRKEMENR RRNQELEEF S KLPFSTVPEA LSFFIARMER DESILEFLQR QSGNKKSYKK KKNNTAEGIS PERKSADAFR KKLIELITAG LTFLEDKIGK EDIYSETRES LQRIYQKLTS NSWSSPVS YD DSNSSQYNFK WEFDDKTYGP YTASQIQAWS NEGYFTDAKH AAFIQLANMD EWMYPNNICF CDVVSLKK
Source	E.coli
Target Names	SPBC83.09c
Protein Names	Recommended name: LIN1-like protein
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.