



# Recombinant Schizosaccharomyces pombe Uncharacterized protein C1604.18c (SPBC1604.18c)

<b>Product Code</b>	CSB-MP530963SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O94730
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGKKEKVNAL EFSVNKGIFT KTRLKSLYSD FTSLFIKNPE GFLANVNTWR EAIETSAWSG KLNSRIILVF DETFESAFSS PALGRPLSLG AVAAYWIIQQG VWLRKQDFLN DCKKGNLVRE DGFSIFSILR WSFQKLGFGQN QASSILSNRS PSGQYVIRKN IEKLACLVHN EAMRRCSSYT SAIYTWDFEQ DTFGSLYWEE GKLNKDEMEC LLNWMCYQKH VLIFDSKIIK FLPNSQIDAQ LASIDKSIDG SVADLIQARA SIAQRSEFLN EELEQLSQVL NQAVKKGEKT IAITYLRRKK ILSKDLERKV SSRLQLDTII SNIDNAVDNK ILLIAMSSGS EALDAILAQM GGTEKVEDVL ENVNDTLARS EEIDATIPTY NPQNIDLEDE AVEKEWQDLV AEEQKVEDIV STLGNVSLKT PSDTFTLTNT DSDKKTSKPE KIQAEELVEQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SPBC1604.18c
<b>Protein Names</b>	Recommended name: Uncharacterized protein C1604.18c
<b>Expression Region</b>	1-449
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.