



Recombinant Drosophila melanogaster Casein kinase II subunit beta' (CkIibeta2)

Product Code	CSB-BP530982DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O96863
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MTDSESSWI HWFCKQRGNE FFCEVDEEYI QDKFNLNFLD SNVKNYKCAL EVILDLNPGS ASEDPAEPEL EASAEKLYGL IHARFILTNR GIELMLDKYN KGEFGTCPRA FCHSQPVLP I GLSDNPGEDM VRIYCPKCND VYIPKASRHS NLDGAFFGTG FPHMFFMEKP DARPKRAKQK FVPRLYGFKI HPTAYRTAAE IQKDVTMTPV GEIDSPSHI
Source	Baculovirus
Target Names	CkIibeta2
Protein Names	Recommended name: Casein kinase II subunit beta' Short name= CK II beta' Short name= DmCKII-beta'
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.