



Recombinant Rabbit Voltage-dependent L-type calcium channel subunit beta-1 (CACNB1), partial

| | |
|--------------------------|---|
| Product Code | CSB-EP004411RB-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P19517 |
| Product Type | Recombinant Protein |
| Immunogen Species | Oryctolagus cuniculus (Rabbit) |
| Purity | ≥85% (SDS-PAGE) |
| Source | E.coli |
| Target Names | CACNB1 |
| Protein Names | Recommended name: Voltage-dependent L-type calcium channel subunit beta-1 Short name= CAB1 Alternative name(s): Calcium channel voltage-dependent subunit beta 1 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Target Details | This protein belongs to the calcium channel beta subunit family. It plays an important role in the calcium channel by modulating G protein inhibition, increasing peak calcium current, controlling the alpha-1 subunit membrane targeting and shifting the voltage dependence of activation and inactivation. Alternative splicing occurs at this locus and three transcript variants encoding three distinct isoforms have been identified. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |