



Recombinant Danio rerio UPF0692 protein C19orf54 homolog (si:dkey-233h2.2, zgc:153119)

Product Code	CSB-YP531777DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B0V3H4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSLENDAAAP PPPPLPPPPP PQPPSLARSE SSKKKLYQAL AEGRTAVEGD YEEASIIIGQ RESSFSKDLQ WLLFNYYVPS LIQDGPQCGL VALWMAAHL L QPPKTLTLLLET LVQTAKDNGY TEQGEMFSAS DMAKLAEDVC GCRVQRLSGG MLGENTAVIL KHLINRQPIL IPYDEDFNHE PCLRNGHKAH WAVASGILLG LREGCVNYKH FPPDITLPWL RLAQSEASVS WPIDDIEEVF VLAKQGKSLR YQLWEFESVA QSNKQLKEMD PQRASDGTRY VLPPGGVQDG LAGQVLLLHT NTEQTKN
Source	Yeast
Target Names	si:dkey-233h2.2
Protein Names	Recommended name: UPF0692 protein C19orf54 homolog
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.