



Recombinant Mouse Thioredoxin-dependent peroxide reductase, mitochondrial (Prdx3)

Product Code	CSB-EP018656MO-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20108
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	TPAVTQHAPYFKGTAVVNGEFKELSLDDFKGKYLVLFFYPLDFTFVCPTEIVAF SDKANEFHDVNCEVVAVSVDSHFSHLAWINTPRKNGGLGHMNITLLSDITKQIS RDYGVLLSAGIALRGLFIIDPNGVVKHLSVNDLPVGRSVEETLRLVKAFQFVET HGEVCPANWTPESPTIKPSPTASKEYFEKVHQ
Source	E.coli
Target Names	Prdx3
Protein Names	Recommended name: Thioredoxin-dependent peroxide reductase, mitochondrial EC= 1.11.1.15 Alternative name(s): Antioxidant protein 1 Short name= AOP-1 PRX III Peroiredoxin-3 Protein MER5
Expression Region	63-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.