



Recombinant Escherichia coli Anaerobic nitric oxide reductase flavorubredoxin (norV)

Product Code	CSB-MP533084ENP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1IUW9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Purity	≥85% (SDS-PAGE)
Sequence	MSIVVKNNIH WVGQRDWEVR DFHGTEYKTL RGSSYNSYLI REEKNVLIDT VDHKFSREFV QNLRNEIDLA DIDYIVINHA EEDHAGALTE LMAQIPDTP I YCTANAIDSI NGH H H H PEWN FNVVKTGDTL DIGNGKQLIF VETPMLHWP D SMMTYLTGDA VLFSNDAFGQ HYCDEHLFND EVDQTELFEQ CQRYANILT PFSRLVTPKI TEILGFNLPV DMIATSHGVV WRDNPTQIVE LYLKWAADYQ EDRTIFYDT MSNTRMMAD AIAQGIAETD PRVAVKIFNV ARSDKNEILT NVFRSKGVLV GTSTMNNVMM PKIAGLVEEM TGLRFRNKRA SAFGSHGWSG GAVDRLSTRL QDAGFEMSLS LKAKWRPDQD ALKLCREHGR EIARQWALAP LPQSTVNTVV KEETSATTTA DLGPRMQCSV CQWIYDPAKG EPMQDVAPGT PWSEVPDNFL CPECSLGKDV FEELASEAK
Source	Mammalian cell
Target Names	norV
Protein Names	Recommended name: Anaerobic nitric oxide reductase flavorubredoxin Short name= FIRd Short name= FlavoRb
Expression Region	1-479
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.