



Recombinant Escherichia coli 4-alpha-L-fucosyltransferase (wecF)

Product Code	CSB-BP534339ENP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1IWA7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Purity	≥85% (SDS-PAGE)
Sequence	MTVLIHVLGS DIPHHNRTVL RFFNDALAAT SEHAREFMVV GKDDGLSDSC PALSQFFPFG KKSLEAVIA KAKANRQQRFFFHGQFNPTL WLALLSGGIK PSQFFWHIWG ADLYELSSGL RYKLFYPLRR LAQKRVGCVF ATRGDLSFFA KTHPKVRGEL LYFPTRMDPS LNTMANDRQR EGKMTILVGN SGDRSNEHVA ALRAVHQQFG DTVKVVVPMG YPPNNEAYIE EVRQAGLELF SEENLQVLSE KLEFDAYLAL LRQCDLGYFI FARQQGIGTL CLLIQAGIPC VLNRENPFWQ DMTEQHLPVL FTTDDLNEDI VREAQRQLAS VDKNTIAFFS PNYLQGWQRA LAIAAGEVA
Source	Baculovirus
Target Names	wecF
Protein Names	Recommended name: 4-alpha-L-fucosyltransferase EC= 2.4.1.- Alternative name(s): TDP-Fuc4NAc:lipid II Fuc4NAc transferase Short name= Fuc4NAc transferase
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.