



# Recombinant Escherichia coli 5-keto-4-deoxy-D-glucarate aldolase (garL)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP535549ENP                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |
| <b>Uniprot No.</b>       | B1IRJ2                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MNNDVFPNKF KAALAAKQVQ IGCWSALSNP ISTEVLGLAG FDWLVL DGEH<br>APNDISTFIP QLMALKGSAS APVVRVPTNE PVIKRLLDI GFYNFLIPFV<br>ETKEKAELAV ASTRYPPEGI RGVSVSHRAN MFGTVADYFA QSNKNITILV<br>QIESQQGVND VDAIAATEGV DGIFVGPSDL AAALGHLGNA SHPDVQKAIQ<br>HIFNRASAHG KPSGILAPVE ADARRYLEWG ATFAVAGSDL GVFRSATQKL<br>ADTFKK                                                                                                    |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Target Names</b>      | garL                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Protein Names</b>     | Recommended name: 5-keto-4-deoxy-D-glucarate aldolase Short name= KDGluc aldolase Short name= KDGlucA EC= 4.1.2.20 Alternative name(s): 2-dehydro-3-deoxy-D-glucarate aldolase 2-keto-3-deoxy-D-glucarate aldolase 5-dehyd                                                                                                                                                                                  |
| <b>Expression Region</b> | 1-256                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |