



Recombinant Escherichia coli High frequency lysogenization protein HfID (hfID)

Product Code	CSB-EP535570ENP-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1IUD5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Purity	>85% (SDS-PAGE)
Sequence	MAKNYYDITL ALAGICQSAR LVQQLAHQGH CDADALHVSL NSIIDMNPSS TLAVFGGSEA NLRVGLETLL GVLNASSRQG LNAELTRYTL SLMVLERKLS SAKGALDTLG NRINGLQRQL EHF DLQSETL MSAMAAIYVD VISPLGPRIQ VTGSPAVLQS PQVQAKVRAT LLAGIRAAVL WHQVGGGRLQ LMFSRNRLTT QAKQILAHLT PEL
Source	E.coli
Target Names	hfID
Protein Names	Recommended name: High frequency lysogenization protein HfID
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.