



# Recombinant Escherichia coli Protein RecA (recA)

<b>Product Code</b>	CSB-BP535577ENP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B1IU0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAIDENKQKA LAAALGQIEK QFGKGSIMRL GEDRSMDVET ISTGSLSLDI ALGAGGLPMG RIVEIYGPES SGKTTTLQV IAAAQREGKT CAFIDAHAL DPIYARKLGV DIDNLLCSQP DTGEQALEIC DALARSGAVD VIVVDSVAAL TPKAEIEGEI GDSHMGLAAR MMSQAMRCLA GNLKQSNTLL IFINQIRMKI GVMFGNPETT TGGNALKFYA SVRLDIRRIG AVKEGENVVG SETRVKVVKN KIAAPFKQAE FQILYEGIN FYGELVDLGV KEKLIEKAGA WYSYKGEKIG QGKANATAWL KDNPETAKEI EKKVRELLLS NPNSTPDFSV DDSEGVAETN EDF
<b>Source</b>	Baculovirus
<b>Target Names</b>	recA
<b>Protein Names</b>	Recommended name: Protein RecA Alternative name(s): Recombinase A
<b>Expression Region</b>	1-353
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.