



Recombinant Escherichia coli Putative transport protein YidE (yidE)

Product Code	CSB-EP535619ENP-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1IYQ9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Purity	>85% (SDS-PAGE)
Sequence	MSDIALTVSILALVAVVGLFIGNVKFRGIGLGIGGVLFGGIIVGHFVSQAGMTLSS DMLH VIQEFGLILFVYTIGIQVGPFFASLRVSGRLRLNLFVAVLIVIGGLVTAILHKLFDIPL P VVLGIFSGAVTNTPALGAGQQILRDLGTPMEMVDQMMSYAMAYPFGICGILF TMWMLRV IFRVNVETEAQQHESSRTNGGALIKTINIRVENPNLHDLAIKDVPILNGDKIICSRL KRE ETLKVPSPDTIIQLGDLLHLVGQPADLHNAQLVIGQEVDTSLSTKGTDLRVERV VVTNEN VLGKRIRDHLHFKERYDVVISRLNRAGVELVASGDISLQFGDILNLVGRPSAIDAV ANVLG NAQQKLQQVQMLPVFIGIGLVLLGSIPVFVPGFPAALKLGLAGGPLIMALILGR IGSIG KLYWFMPPSANLALRELGIVLFLSVVGLKSGGDFVNTLVNGEGLSWIGYGALIT AVPLIT VGILARMLAKMNYLTMCGMLAGSMTDPPALAFANNLHPTSGAAALSYATVYPL VMFLRII TPQLLAVLFWSIG
Source	E.coli
Target Names	yidE
Protein Names	Recommended name: Putative transport protein YidE
Expression Region	1-553
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.