



# Recombinant Escherichia coli Ribulokinase (araB)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP536799ENP  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | B11RB5   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | MAIAIGLDFG SDSVRALAVD CASGEEIATS VEWYPRWQKG QFCDAPNNQF<br>RHHPRDYIES MEAALKTVLA ELSVEQRAAV VGIGVDSTGS TPAPIDADGN<br>VLALRPEFAE NPNAMFVLWK DHTAVEEAE ITRLCHAPGN VDYSRYIGGI<br>YSSEFWAKI LHVTRQDSAV AQAASWIEL CDWVPALLSG TTRPQDIRRG<br>RCSAGHKSLW HESWGGLPPA SFFDELDPIL NRHLPSPFLT DTWTADIPVG<br>TLCPEWAQRL GLPESVVISG GAFDCHMGAV GAGAQPNALV KVIGTSTCDI<br>LIADKQSVGE RAVKGICGQV DGSVVPFGFIG LEAGQSAFGD IYAWFGRVLS<br>WPLEQLAAQH PELKAQINAS QKQLLPALTE AWAKNPSLDH LPVVLDFWNG<br>RRSPNANQRL KGVITDLNLA TDAPLLFGGL IAATAFGARA IMECFTDQGI<br>AVNNVMALGG IARKNQVIMQ ACCDVLNRPL QIVASDQCCA LGAAIFAAVA<br>AKVHADIPSA QKMASAVEK TLQPRSEQAQ RFEQLYRRYQ<br>QWAMSAEQHY LPTSAPAQAA QAVATL |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | araB   |
| <b>Protein Names</b>     | Recommended name: Ribulokinase EC= 2.7.1.16  |
| <b>Expression Region</b> | 1-566  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |