



Recombinant Escherichia coli Fructose-1,6-bisphosphatase class 1 (fbp)

Product Code	CSB-EP536817ENP-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1ISX3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Purity	>85% (SDS-PAGE)
Sequence	MKTLGEFIVE KQHEFSHATG ELTALLSAIK LGAKIIHRDI NKAGLVDILG ASGAENVQGE VQQKLDLFAN EKLKAALKAR DIVAGIASEE EDEIVVFEGC EHAKYVVLMD PLDGSSNIDV NVSVGTIFSI YRRVTPVGTP VTEEDFLQPG NKQVAAGYVV YGSSTMLVYT TGCGVHAFTY DPSLGVFCLC QERMRFPEKG KTYSINEGNY IKFPNGVKKY IKFCQEEDKS TNRPYTSRYI GSLVADFHRN LLKGGIYLYP STASHPDGKL RLLYECNPMA FLAEQAGGKA SDGKERILDI IPETLHQRRS FFGNDHMVE DVERFIREFP DA
Source	E.coli
Target Names	fbp
Protein Names	Recommended name: Fructose-1,6-bisphosphatase class 1 Short name=FBPase class 1 EC= 3.1.3.11 Alternative name(s): D-fructose-1,6-bisphosphate 1-phosphohydrolase class 1
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.