



# Recombinant Escherichia coli Sulfite reductase [NADPH] hemoprotein beta-component (cysl)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP536826ENP   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B1IU78  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSEKHPGPLV VEGKLTDAER MKLESNYLRG TIAEDLNDGL TGGFKGDNFL LIRFHGMYYQ DDRDIRAERA EQKLEPRHAM LLRCRLPGGV ITTKQWQAID KFAAGENTIYG SIRLTNRQTF QFHGILKKNV KPVHQLHSV GLDALATAND MNRNVLCTSN PYESQLHAEA YEWAKKISEH LLPRTRAYAE IWLDQEKVAT TDEEPILGQT YLPRKFKTTV VIPPQNDIDL HANDMNFVAI AENGKLVGFN LLVGGGLSIE HGNKKTYART ASEFGYLPLE HTLAVAEAVV TTQRDWGNRT DRKNAKTKYT LERVGVEFK AEVERRAGIK FEPIRPYEFT GRGDRIGWVK GIDDNWHLTL FIENGRILDY PGRPLKTGLL EIAKIHKGDF RITANQNLI AGVPESEKAK IEKIAKESGL MNAVTPQREN SMACVSFPTC PLAMAEAEERF LPSFIDNIDN LMAKHGVSDE HIVMRVTGCP NGCGRAMLAE VGLVGKAPGR YNLHLSGNRI GTRIPRMYKE NITEPEILAS LDELIGRWAK EREAGEGFGD FTVRAGIIRP VLDPARDLWD |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | cysl  |
| <b>Protein Names</b>     | Recommended name: Sulfite reductase [NADPH] hemoprotein beta-component<br>Short name= SiR-HP Short name= SiRHP EC= 1.8.1.2  |
| <b>Expression Region</b> | 1-570   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.