



# Recombinant Escherichia coli Allantoinase (allB)

<b>Product Code</b>	CSB-EP536879ENP-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B11Z89
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSFDLIKNG TVILENEARV VDIADVGGKI AAIGQDLGDA KEVMDASGLV VSPGMVDAHT HISEPGRSHW EGYETGTRAA AKGGITTMIE MPLNQLPATV DRASIELKFD AAKGKLTIDA AQLGGLVSYN IDRLHELDEV GVVGFKCFVA TCGDRGIDND FRDVNDWQFF KGAQKLGELG QPVLVHCENA LICDALGEEA KSEGRVTAHD YVASRPVFTE VEAIRRVLYL AKVAGCRLHI CHISSPEGVE EVTRARQEGQ DVTCESEPHY FVLDTDQFEE IGTLAKCSPP IRDLNQKGM WEKLFNGEID CLVSDHSPCP PEMKAGNIME AWGGIAGLQN CMDVMFDEAV QKRGMSLPMF GKLMATNAAD IFGLQQKGRI APGKDADFVF IQPNSSYVLT NDDLEYRHKV SPYVVRTIGA RITKTILRGD VIYDIEQGFP VAPKGQFILK HQQ
<b>Source</b>	E.coli
<b>Target Names</b>	allB
<b>Protein Names</b>	Recommended name: Allantoinase EC= 3.5.2.5 Alternative name(s): Allantoin-utilizing enzyme
<b>Expression Region</b>	1-453
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.