



Recombinant Escherichia coli Gamma-aminobutyraldehyde dehydrogenase (prp)

Product Code	CSB-EP538052ENP-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1IS15
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Purity	≥85% (SDS-PAGE)
Sequence	<pre>MQHKLLINGE LVSgegekqp VYNPATGDVL LEIAEASAEQ VDAAVRAADA AFAEWGQ TTP KVRAECLLKL ADVIEENGQV FAELESRNCG KPLHS AFNDE IPAIVDVFRF FAGAARCLNG LAAGEYLEGH TSMIRRDPLG VVASIAPWNY PLMMAAWKLA PALAAGNCVV LKPSEITPLT ALKLAELAKD IFPAGVINIL FGRGKTVGDP LTGHPKVRMV SLTGS IATGE HIISHTASSI KRTHMELGGK APVIVFDDAD IEAVVEGVRT FGYYNAGQDC TAACRIYAQK GIYDTLVEKL GAAVATLKSG APDDESTELG PLSSLAHLER VSKAVEEAKA TGHKIVITGG EK RKGN GYYY APTLLAGALQ DDAIVQKEVF GPVVS VTLFD NEEQV VNWAN DSQYGLASSV WTKDVGRAHR VSARLQYGCT WVNTHFMLVS EMPHGGQKLS GYGKDMSLYG LEDYTVVRHV MVKH</pre>
Source	E.coli
Target Names	prp
Protein Names	Recommended name: Gamma-aminobutyraldehyde dehydrogenase EC= 1.2.1.19 Alternative name(s): 1-pyrroline dehydrogenase 4-aminobutanal dehydrogenase Short name= ABALDH
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.