



# Recombinant Escherichia coli L-fucose isomerase (fucI)

<b>Product Code</b>	CSB-MP538064ENP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B1IU39
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKKISLPKIG IRPVIDGRRM GVRESLEEQT MNMAKATAAL LTEKLRHACG AAVECVISDT CIAGMAEAAA CEEKFSSQNV GLTITVTPCW CYGSETIDMD PTRPKAIWGF NGTERPGAVY LAAALAAHSQ KGIPAFSIYG HDVQDADDT IPADVEEKLL RFARAGLAVA SMKGKSYLSL GGVSMMGIAGS IVDHNFESW LGMKVQAVDM TELRRRIDQQ IYDEAELEMA LAWADKNFRY GEDENNKQYQ RNAEQSRAVL RESLLMAMCI RDMMQGNLADIGRVEESL GYNAIAAGFQ GQRHWTDQYP NGDTAEAILN SSFDWNGVRE PFVVATENDS LNGVAMLMGH QLTGTAQVFA DVRTYWSPEA IERVTGHKLD GLAEHGIHL INSGSAALDG SCKQRDSEGN PTMKPHWEIS QQEADACLAA TEWCPAIHEY FRGGGYSSRF LTEGGVPFTM TRVNIKGLG PVLQIAEGWS VELPKDVHDI LNKRTNSTWP TTWFAPRLTG KGPFTDVYSV MANWGANHGV LTIGHVGADF ITLASMLRIP VCMHNVEETK VYRPSAWAAH GMDIEGQDYR ACQNYGPLYK R
<b>Source</b>	Mammalian cell
<b>Target Names</b>	fucI
<b>Protein Names</b>	Recommended name: L-fucose isomerase EC= 5.3.1.25 Alternative name(s): 6-deoxy-L-galactose isomerase FucIase
<b>Expression Region</b>	1-591
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.