



# Recombinant Escherichia coli tRNA modification GTPase MnmE (mnmE)

<b>Product Code</b>	CSB-BP538103ENP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B1IX32
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSDNDTIVAQ ATPPGRGGVG ILRISGFKAR EVAETVLGKL PKPRYADYLP FKDADGSVLD QGIALWFPGP NSFTGEDVLE LQGHGGPVIL DLLLKRILTI PGLRIARPGE FSERAF LNDK LD LAQAE AIA DLIDASSEQA ARSALNSLQG AFSARVNHLV EALTHLRIYV EAAIDFPDEE IDFLSDGKIE AQLNDVIADL DAVRAEARQG SLLREGMKVV IAGRPNAGKS SLLNALAGRE AAIVTDIAGT TRDVLREHIIH IDGMPLHIID TAGLREASDE VERIGIERAW QEIEQADRVL FMVDGTTTDA VDP AEIWPEF IARLPAKLPI TVVRNKADIT GETLGMSEVN GHALIRLSAR TGEGVDVLRN HLKQSMGFD T NMEGGFLARR RHLQALEQAA EHLQQGKAQL LGAWAGELLA EELRLAQQNL SEITGEFTSD DLLGRIFSSF CIGK
<b>Source</b>	Baculovirus
<b>Target Names</b>	mnmE
<b>Protein Names</b>	Recommended name: tRNA modification GTPase MnmE EC= 3.6.-.-
<b>Expression Region</b>	1-454
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.