



# Recombinant Escherichia coli S-(hydroxymethyl)glutathione dehydrogenase

<b>Product Code</b>	CSB-MP538132ENP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B1J085
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKSRAAVAFA PGKPLEIVEI DVAPPKKGEV LIKVHTHTGVC HTDAFTLSGD DPEGVFPVVL GHEGAGVVVE VGEGVTSVKP GDHVIPLYTA ECGECEFCRS GKTNLCVAVR ETQGKGLMPD GTTRFSYNGQ PLYHYMGCST FSEYTVVAEV SLAKINPEAN HEHVCLLGCG VTTGIGAVHN TAKVQPGDSV AVFGLGAIGL AVVQGARQAK AGRIIAIDTN PKKFDLARRF GATDCINPND YDKPIKDVL DINKWGIDHT FECIGNVNVM RAALESAHRG WGQSVIIGVA GSGQEISTRP FQLVTGRVWK GSAFGGVKGR SQLPGMVEDA MKGDIDLEPF VTHTMSLDEI NDAFDLMHEG KSIRTVIRY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	frmA
<b>Protein Names</b>	Recommended name: S-(hydroxymethyl)glutathione dehydrogenase EC= 1.1.1.284 Alternative name(s): Alcohol dehydrogenase class-3 EC= 1.1.1.1 Alcohol dehydrogenase class-III Glutathione-dependent formaldehyde dehydrogenase Sh
<b>Expression Region</b>	1-369
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.