



Recombinant Human Protein-L-isoaspartate (D-aspartate) O-methyltransferase (PCMT1)

Product Code	CSB-EP017618HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P22061
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AWKSGGASH SELIHNLRKN GIIKTDKVEE VMLATDRSHY AKCNPYMDSP QSIGFQATIS APHMHAYALE LLFDQLHEGA KALDVGSGSG ILTACFARMV GCTGKVIDID HIKELVDDSV NNVRKDDPTL LSSGRVQLVV GDGRMGYAE APYDAIHVGA AAPVVPQALI DQLKPGGRLI LPVGPAGGNQ MLEQYDKLQD GSIKMKPLMG VIYVPLTDKE KQWSRWK
Source	E.coli
Target Names	PCMT1
Protein Names	Recommended name: Protein-L-isoaspartate(D-aspartate) O-methyltransferase Short name= PIMT EC= 2.1.1.77 Alternative name(s): L-isoaspartyl protein carboxyl methyltransferase Protein L-isoaspartyl/D-aspartyl methyltransferase
Expression Region	2-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.