



Recombinant Escherichia coli Threonine--tRNA ligase (thrS)

Product Code	CSB-MP539298ENP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1IPK9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Purity	>85% (SDS-PAGE)
Sequence	MPVITLPDGS QRHYDHAVSP MDVALDIGPG LAKACIAGRV NGELVDACDL IENDAQLSII TAKDEDGLEI IRHSCAHLG HAIKQLWPHT KMAIGPVIDN GFYYDVDLDR TLTQEDVEAL EKRMHELAEK NYDVIKKKVS WHEARETFAN RGESYKVSIL DENIAHDDKP GLYFHEEYVD MCRGPHVPMN RFCHHFKLMK TAGAYWRGDS NNKMLQRIYG TAWADKKALN AYLQRLEEEA KRDRHKIGKQ LDLYHMQEEA PGMVFWHNDG WTIFRELEV FVRSKLKEYQY QEVKGPMMMD RVLWEKTGHW DNYKDAMFTT SSENREYCIK PMNCPGHVQI FNQGLKSYRD LPLRMAEFGS CHRNEPSGSL HGLMRVRGFT QDDAHIFCTE EQIRDEVNGC IRLVYDMYST FGFEKIVVKL STRPEKRIGS DEMWDRAEAD LAVALEENNI PFEYQLGEGA FYGPKIEFTL YDCLDRAWQC GTVQLDFSLP SRLSASYVGE DNERKVPVMI HRAILGSMER FIGILTEEFA GFFPTWLAPV QVIMNITDS QSDYVNETQ KLSNAGIRVK ADLRNEKIGF KIREHTLRRV PYMLVCGDKE VESGKVAVRT RRGKDLGSMD VNEVIEKLQQ EIRSRSLKQL EE
Source	Mammalian cell
Target Names	thrS
Protein Names	Recommended name: Threonine--tRNA ligase EC= 6.1.1.3 Alternative name(s): Threonyl-tRNA synthetase Short name= ThrRS
Expression Region	1-642
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.