



# Recombinant Escherichia coli Cytochrome c-552 (nrfA)

<b>Product Code</b>	CSB-YP539349ENP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B1IU11
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>EQTA APAKPVTVEA KNETFAPQHP DQYLSWKATS EQSERVDALA  EDPRLVILWA GYPFSRDYNK PRGHAFVTD VRETLRTGAP KNAEDGPLPM  ACWSCKSPDV ARLIQKDGED GYFHGKWARG GPEIVNNLGC ADCHNTASPE  FAKGKPELTL SRPYAARAME AIGKPFEEKAG RFDQQSMVCG QCHVEYYFDG  KNKAVKFPWD DGMKVENMEQ YYDKIAFSDW TNSLSKTPML  KAQHPEYETW TAGIHGKNNV TCIDCHMPKV QNAEGKLYTD HKIGNPFDNF  AQTCANCHTQ DKAALQKVVA ERKQSINDLK IKVEDQLVHA HFEAKAALDA  GATEAEMKPI QDDIRHAQWR WDLAIASHGI HMHAPEEGLR MLGTAMDKAA  DARTKLARLL ATKGITHEIQ IPDISTKEKA QQAIGLNMEQ IKAEEKQDFIK  TVIPQWEEQA RKNGLLSQ</p>
<b>Source</b>	Yeast
<b>Target Names</b>	nrfA
<b>Protein Names</b>	Recommended name: Cytochrome c-552 EC= 1.7.2.2 Alternative name(s): Ammonia-forming cytochrome c nitrite reductase Short name= Cytochrome c nitrite reductase
<b>Expression Region</b>	27-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.