



# Recombinant Escherichia coli Fatty acid oxidation complex subunit alpha (fadJ)

<b>Product Code</b>	CSB-YP539382ENP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B1IXA5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEMASAFTLN VRLDNIAIIT IDVPDEKMNT LKAEFASQVR AIIKQLRENK ELRGVVFISA KPDNFIAGAD INMIGNCKTA QEAEALARQG QQLMAEIHAI PIPVIAAIHG ACLGGGLELA LACHGRVCTD DPKTVLGLPE VQLGLLPGSG GTQRLPRLIG VSTALEMILT GKQLRAKQAL KLGLVDDVVP HSILLEAAVE LVKQDRPSSR PLPVRERILA GPLGRALLFK MVGKKTEHKT QGNYPATERI LEVVETGLAQ GTSSGYDAEA RAFGELAMTP QSQUALRSIFF ASTDVKKDPG SDAPPAPLNS VGILGGGLMG GGIAVVTACK AGLPVRIKDI NPQGINHALK YSWDQLEGKV RRRHLKASER DKQLALISGT TDYRGFAHRD LIIEAVFENL ELKQQMVAEV EQNCAAHTIF ASNTSSLPIG DIAAHATRPE QVIGLHFFSP VEKMPLVEII PHAGTSAQTI ATTVKLAKKQ GKTPIVVRDK AGFYVNRILA PYINEAIRML TEGERVEHID AALVKFGFPV GPIQLLDEVG IDTGTKIIPV LEAAYGERFS APANVVSSIL NDDRKGRKNG RGFYLYGQKG RSKKKQVDPA IYPLIGAQQG GRLSAPQVAE RCVMLMLNEA VRCVDEQVIR SVRDGDIGAV FGIGFPPFLG GPFYIDSLG AGEVVAIMQR LATQYGSRFT PCERLVEMSK RGESFWKTTA TDLQ
<b>Source</b>	Yeast
<b>Target Names</b>	fadJ
<b>Protein Names</b>	Recommended name: Fatty acid oxidation complex subunit alpha Including the following 2 domains: Enoyl-CoA hydratase/3-hydroxybutyryl-CoA epimerase EC= 4.2.1.17 EC= 5.1.2.3 3-hydroxyacyl-CoA dehydrogenase EC= 1.1.1.35
<b>Expression Region</b>	1-714
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.