



Recombinant Escherichia coli tRNA 5-methylaminomethyl-2-thiouridine biosynthesis bifunctional protein MnmC (mnmC)

Product Code	CSB-BP539383ENP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1IXC1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Purity	≥85% (SDS-PAGE)
Sequence	MKHYSIQPAN LEFNAEGTPV SRDFDDVYFS NDNGLEETRY VFLGGNQLEV RFPEHPHPLF VVAESGFGTG LNFLTLLWQAF DQFREAHPQA RLQRLHFISF EKFPPLTRADL ALAHQHWPEL APWAEQLQAQ WPMPLPGCHR LLLDEGRVTL DLWFGDINEL TSQLDDSLNQ KVDAWFLDGF APAKNPDMWT QNLFNAMARL ARPGGTLATF TSAGFVRRGL QEAGFTMQKR KGFGGRKREML CGVMEQTLPL PCSAPWFNRT GSSKREAAII GGGIASALLS LALLRRGWQV TLYCADEAPA LGASGNRQGA LYPLLSKHDE ALNRFFSNAF TFARRFYDQL PVKFDHDWCG VTQLGWDEKS QHKIAQMLSM DLPAELAVAV EANAVEQITG VATNCSGITY PQGGWLCPAE LTRNVLELAQ QQGLQIYYQY QLQNLSRKDD CWLLNFAGDQ QATHSVVLA NGHQISRFSQ TSTLPVYSVA GQVSHIPTTP ELAELKQVLC YDGYLTPQNP ANQHHICIGAS YHRGSEDYAY SEDDQQQNRQ RLIDCFPQAQ WAKEVDVSDK EARCGVRCAT RDHLPVGNV PDYEATLVEY ASLAEQKDEA VSAPVFDDL F MFAALGSRGL CSAPLCAEIL AAQMSDEPIP MDASTLAALN PNRLWVRKLL KGKAVKAG
Source	Baculovirus
Target Names	mnmC
Protein Names	Recommended name: tRNA 5-methylaminomethyl-2-thiouridine biosynthesis bifunctional protein MnmC Short name= tRNA mnm(5)s(2)U biosynthesis bifunctional protein Including the following 2 domains: tRNA (mnm(5)s(2)U34)-methyltransferase
Expression Region	1-668
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.