



Recombinant Escherichia coli Chorismate synthase (aroC)

Product Code	CSB-YP540365ENT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1X9K2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12 / DH10B)
Purity	>85% (SDS-PAGE)
Sequence	MAGNTIGQLF RVTTFGESHG LALGCIVDGV PPGIPLTEAD LQHDLDRRRP GTSRYTTQRR EPDQVKILSG VFEGVTTGTS IGLLIENTDQ RSQDYSIAKD VFRPGHADYT YEQKYGLRDY RGGGRSSARE TAMRVAAGAI AKKYLAEKFG IEIRGCLTQM GDIPLDIKDW SQVEQNPFFC PDPDKIDALD ELMRALKKEG DSIGAKVTVV ASGVPAGLGE PVFDRDLADI AHALMSINAV KGVEIGDGF VVALRGSQNR DEITKDFQSQ NHAGGILGGI SSGQQIAHM ALKPTSSITV PGRTINRFGE EVEMITKGRH DPCV GIRAVP IAEAMLAI VL MDHLLRQRAQ NADVKT DIPR W
Source	Yeast
Target Names	aroC
Protein Names	Recommended name: Chorismate synthase EC= 4.2.3.5 Alternative name(s): 5-enolpyruvylshikimate-3-phosphate phospholyase
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.