



# Recombinant Escherichia coli Putative ribosome biogenesis GTPase RsgA (rsgA)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP540412ENT-B  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | B1XDR5   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli (strain K12 / DH10B)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MSKNKLSKGQ QRRVNANHQR RLKTSKEKPD YDDNLFGEPE EGVISRFGM<br>HADVESADGD VHR CNIRRTI RSLVTGDRVV WRPGKPA AEG VNVKGIVEAV<br>HERTSVLTRP DFYDGVKPIA ANIDQIVIVS AILPELSLNI IDRYLVACET<br>LQIEPIIVLN KIDLLDDEGM AFVNEQMDIY RNIGYRVLMV SSHTQDGLKP<br>LEEAL TGRIS IFAGQSGVGK SLLNALLGL QKEILTNDIS DNSGLGQHTT<br>TAARLYHFPH GGDVIDSPGV REFGLWHLEP EQITQGFVEF HDYLG LCKYR<br>DCKHDTDPGC AIREAVEEGK IAETRFENYH RILESMAQVK TRKNFSDTDD |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | rsgA   |
| <b>Protein Names</b>     | Recommended name: Putative ribosome biogenesis GTPase RsgA EC= 3.6.1.-   |
| <b>Expression Region</b> | 1-350  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.          |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |