



Recombinant Escherichia coli Probable endonuclease 4 (nfo)

Product Code	CSB-EP540724ENX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1LKQ6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain SMS-3-5 / SECEC)
Purity	>85% (SDS-PAGE)
Sequence	MKYIGAHVSA AGGLANAAIR AAEIDATAFA LFTKNQRQWR AAPLTTQTID EFKAACEKYH YTSAQILPHD SYLINLGHVPV TEALEKSRDA FIDEMQRCEQ LGLSLLNFHP GSHLMQISEE DCLARIAESI NIALDKTHGV TAVIENTAGQ GSNLGFKFEH LAAIIDGVED KSRVGVICIDT CHAFAAGYDL RTSAECEKTF ADFARIVGFK YLRGMHLNDA KSTFGSRVDR HHSLGEGNIG HDAFRWIMQD DRFDGIPLIL ETINPDIWAE EIAWLKAQQT EKAVA
Source	E.coli
Target Names	nfo
Protein Names	Recommended name: Probable endonuclease 4 EC= 3.1.21.2 Alternative name(s): Endodeoxyribonuclease IV Endonuclease IV
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.