



Recombinant Escherichia coli UPF0271 protein YbgL (ybgL)

Product Code	CSB-MP540729ENX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1LLF0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain SMS-3-5 / SECEC)
Purity	>85% (SDS-PAGE)
Sequence	MKIDLNADLG EGYASDAELL TLVSSANIAC GFHAGDAQTM QACVREAIKN GVAIGAHPSF PDRENFGRSA MQLPPETAFA QTLYQIGALA AIARAQGGVM RHVKPHGMLY NQAAKEAQLA DAIARAVYAC DPALILVGLA GSELIRAGKQ YGLTTREEVF ADRGYQADGS LVPRSQPGAL IEDEEQALAQ TLEMVQHGRV KSITGEWATV TAQTVCLHGD GEHALAFARR LRATFAEKGI VVAA
Source	Mammalian cell
Target Names	ybgL
Protein Names	Recommended name: UPF0271 protein YbgL
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.