



Recombinant Mouse Mast cell carboxypeptidase A (Cpa3)

Product Code	CSB-MP005877MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15089
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	I AGRHSYAKYN DWDKIVSWTE KMLEKHPPEMV SRIKIGSTVE DNPLYVLKIG KKDGERKAIF MDCGIHAREW ISPAFCQWFV YQATKSYGKN KIMTKLLDRM NFYVLPVFNV DGYIWSWTQD RMWRKNRSRN QNSTCIGTDL NRNFDVSWDS SPNTNKPCLN VYRGPAPASE KETKAVTNFI RSHLNSIKAY ITFHSYSQML LIPYG YTFKL PPNHQDLLKV ARIATDALST RYETRYIYGP IASTIYKTSG SSLDWVYDLG IKHTFAFELR DKGKSGFLLP ESRIKPTCKE TMLSVKFAK YILKNTS
Source	Mammalian cell
Target Names	Cpa3
Protein Names	Recommended name: Mast cell carboxypeptidase A Short name= MC-CPA EC= 3.4.17.1 Alternative name(s): Carboxypeptidase A3
Expression Region	110-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Three different forms of human pancreatic procarboxypeptidase A have been isolated. This gene encodes a form which is obtained as a binary complex of a procarboxypeptidase A with proproteinase E and functions as a secretory granule metalloexopeptidase.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.