



Recombinant Escherichia coli Glycerol-3-phosphate dehydrogenase [NAD (P)+]

Product Code	CSB-MP541015ENT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1X942
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12 / DH10B)
Purity	>85% (SDS-PAGE)
Sequence	MNQRNASMTV IGAGSYGTAL AITLARNGHE VVLWGHDPPEH IATLERDRCN AAFLPDVFPF DTLHLESDLA TALAASRNIL VVPSHFVGE VLRQIKPLMR PDARLVWATK GLEAETGRLL QDVAREALGD QIPLAVISGP TFAKELAAGL PTAISLASTD QTFADDLQQL LHCCKSFRVY SNPDFIGVQL GGAVKNVIAI GAGMSDGIGF GANARTALIT RGLAEMSRLG AALGADPATF MGMAGLGDLV LTCTDNQSRN RRFGMMLGQG MDVQSAQEKI GQVVEGYRNT KEVRELAHRF GVEMPITEEI YQVLYCGKNA REAALTLLGR ARKDERSSH
Source	Mammalian cell
Target Names	gpsA
Protein Names	Recommended name: Glycerol-3-phosphate dehydrogenase [NAD(P)+] EC=1.1.1.94 Alternative name(s): NAD(P)H-dependent glycerol-3-phosphate dehydrogenase
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.