



Recombinant Callicebus moloch High mobility group protein B1 (HMGB1)

Product Code	CSB-BP541463CYE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1MTB0
Product Type	Recombinant Protein
Immunogen Species	Plecturocebus moloch (Dusky titi monkey) (Callicebus moloch)
Purity	≥85% (SDS-PAGE)
Sequence	GKGDPPKPR GKMSYAFFV QTCREEHKKK HPDASVNFSE FSKKCSERWK TMSAKEKGKF EDMAKADKAR YEREMKTYIP PKGETKKKFK DPNAPKRPPS AFFLFCSEYR PKIKGEHPGL SIGDVAKKLG EMWNNTAADD KQPYEKKA LKEKYEKDIA AYRAKGPDA AKKGVVKA EK SKKKKEEEED EEEDEDEEEE EDEDEDEEEE DDDDE
Source	Baculovirus
Target Names	HMGB1
Protein Names	Recommended name: High mobility group protein B1 Alternative name(s): High mobility group protein 1 Short name= HMG-1
Expression Region	2-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.