



Recombinant Escherichia coli Iron-binding protein IscA (iscA)

| | |
|--------------------------|---|
| Product Code | CSB-EP541693ENT-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B1XB03 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli (strain K12 / DH10B) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MSITLSDSAA ARVNTFLANR GKGFGRLRLGV RTSGCSGMAY VLEFVDEPTP EDIVFEDKGV KVVVDGKSLQ FLDGTQLDFV KEGLNEGFKF TNPNVKDECG CGESFHV |
| Source | E.coli |
| Target Names | iscA |
| Protein Names | Recommended name: Iron-binding protein IscA Alternative name(s): Iron-sulfur cluster assembly protein |
| Expression Region | 1-107 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |