



# Recombinant Escherichia coli Molybdenum cofactor biosynthesis protein A (moaA)

<b>Product Code</b>	CSB-EP542068ENX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B1LM72
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain SMS-3-5 / SECEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASQLTDAFA RKFYYLRLSI TDVCNFRCTY CLPDGYKPSG VTNKGFLTVD EIRRVTAF A SLGTEKVRLT GGEP SLRRDF TDIIAAVREN DAIRQIAVTT NGYRLERDVA NWRDAGLTGI NVSVDSL DAR QFHAI TGQDK FNQVMAGIDA AFEAGFEKVK VNTVLMRDVN HHQLDTFLNW IQHRPIQLRF IELMETGEGS ELFRKHHISG QVLRDELLRR GWIHQ LRQRS DGPAQV FCHP DYAGEIGLIM PYEKDFCATC NRLRVSSIGK LHLCLFGEGG VNLRDLLEDD TQQQALEARI SAALREKKQT HFLHQ NNTGI TQNLSYIGG
<b>Source</b>	E.coli
<b>Target Names</b>	moaA
<b>Protein Names</b>	Recommended name: Molybdenum cofactor biosynthesis protein A
<b>Expression Region</b>	1-329
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.